






Expanding Horizons for the Visually Impaired

Frame Tech Sheet				
 <p>Ocutech Zyl Standard Frame For use with VES II</p>	Eye	DBL	Temples	Max Seg. Ht.
	50	26	140 & 150 mm	6mm
	54	26	140 & 150 mm	6mm
	<i>Standard Colors- Grey, Crystal, Demi Amber, Tan</i>			
	Eye	DBL	Temples	Max Seg. Ht.
46	20	135mm	NR	
<i>Pediatric Colors- Grey, Crystal, Demi Amber, Tan</i>				
 <p>Ocutech Zyl Unisex Frame For use with VES II</p>	Eye	DBL	Temples	Max Seg. Ht.
	48	23.5	140mm	NR
	<i>Unisex Colors- Grey, Crystal, Demi Amber, Demi Blonde</i>			
 <p>Ocutech "K" Metal Frame For use with VES AF, K, Sport, Mini</p>	Eye	DBL	Temples	Max Seg. Ht.
	51	18	135mm	14mm
	53	18	140mm	16mm
	55	18	140mm	19mm
	57	18	145mm	21mm
	59	18	145mm	23mm
<i>"K" Metal Colors- Gold, Natural, Bronze</i>				
 <p>Ocutech Unisex Metal Frame For use with VES AF, K, Sport, Mini</p>	Eye	DBL	Temples	Max Seg. Ht.
	49	18	135mm	9mm
	51	18	140mm	12mm
	<i>Unisex Metal Colors- Gold, Natural, Bronze</i>			
 <p>Ocutech Sleek Metal Frame For use with VES K, Sport, Mini</p>	Eye	DBL	Temples	Max Seg. Ht.
	50	18	140mm	NR
	53	18	145mm	NR
	<i>Sleek Metal Colors- Natural, Gun Metal</i>			